

In the presence of 5% blood serum and a contact time of 5 minutes, Contec® Citric Acid Disinfectant is effective against: **GERMICIDAL:** *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Salmonella enterica*, *Listeria monocytogenes*, *Escherichia coli* O157:H7, *Staphylococcus aureus* - MRSA, *Staphylococcus epidermidis* - MRSE, *Enterococcus faecalis* Vancomycin Resistant (VRE) **VIROCIDAL*:** Respiratory Syncytial Virus, Rotavirus, Human Immunodeficiency Virus Type 1 (Strain HTLV-III₁), Herpes Simplex Type 1 Virus, Herpes Simplex Type 2 Virus (Strain G), Vaccinia Virus, Influenza A Virus (Strain Hong Kong), Pandemic 2009 H1N1 Influenza A Virus (formerly called swine flu), Adenovirus Type 2, Feline Calicivirus (Surrogate for Norovirus) **FUNGICIDAL:** *Trichophyton mentagrophytes*

In the presence of 5% blood serum and a contact time of 10 minutes, Contec® Citric Acid Disinfectant is effective against: **GERMICIDAL:** *Acinetobacter baumannii* ATCC 19606 **VIROCIDAL:** Canine Parvovirus (CPV), Hepatitis B Virus (HBV), Poliovirus Type 1 **TUBERCULOCIDAL:** Effective against *Mycobacterium bovis* (BCG) in 5 minutes at 20°C/68°F on previously cleaned hard non-porous inanimate surfaces. **SANITIZER:** *Staphylococcus aureus* and *Klebsiella pneumoniae* in 1 minute.

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Product is ready to use. **DO NOT DILUTE.** Preclean areas with gross filth and heavy soil prior to disinfection.

TO DEODORIZE: Apply to surfaces as needed. **TO CLEAN:** Apply Contec Citric Acid Disinfectant undiluted to hard, non-porous surfaces by wiping with a cloth or towel until clean (Spot check to test for compatibility). **TO SANITIZE:** Apply to the surface. Let stand for 1 minute. Wipe or let air dry. No rinsing or wiping is required. **TO DISINFECT:** Apply Contec Citric Acid Disinfectant undiluted to hard, non-porous surfaces with a cloth, mop or brush, or by immersing equipment (item) to be disinfected. Wet all surfaces thoroughly. The surface must remain wet with Contec Citric Acid Disinfectant for 10 minutes. Wipe or let air dry. No rinsing or wiping is required, except on direct food contact surfaces which require a potable water rinse after treatment. **Toys:** Use only on hard non-porous surfaces. Rinse toys with potable water. **FUNGICIDAL DIRECTIONS:** This product is effective against *Trichophyton mentagrophytes* (athlete's foot fungus and a cause of ringworm) on inanimate surfaces in locker rooms, dressing rooms, exercise facilities, and shower and bath areas. Follow disinfection directions.

Areas of Use: Athletic facilities, nursing homes, health care facilities, bathrooms, medical offices, institutions, dialysis units, kitchens, daycare centers, grooming facilities, residential, prisons, schools, laboratories or other small animal facilities, restoration, decontamination and remediation sites (fire, water, sewage, trauma scene decontamination and remediation).

Types of surfaces/use sites: Baked enamel, chrome, glass, metal, plastic, stainless steel, and vinyl. For a list of more surfaces/use sites, visit our website www.contecinc.com.

DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION: This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

***CONTEC CITRIC ACID DISINFECTANT KILLS HIV ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus). PERSONAL PROTECTION: When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask and eye coverings.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

CONTACT TIME: Do Not Dilute. Leave precleaned surfaces wet for 5 minutes using mop, cloth or sponge to mitigate HIV-1 virus. Use a 10 minute contact time to kill all bacteria, viruses* and fungi listed on the label.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other fluids should be autoclaved and disposed of according to Federal, State and Local regulations for infectious waste disposal.

VETERINARY PRACTICE/ANIMAL CARE/ANIMAL LABORATORY/ZOOS/PET SHOP/KENNELS DISINFECTION DIRECTIONS: For cleaning and disinfecting hard nonporous surfaces: equipment used for animal food or water, utensils, instruments, cages, kennels, stables, catteries, etc. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings and manure from floors, walls and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Saturate surfaces with Contec Citric Acid Disinfectant for a period of 10 minutes. Wipe or allow to air dry. Thoroughly scrub all treated surfaces that come into contact with food, including equipment used for feeding or watering, with soap or detergent and rinse with potable water before reuse. Do not house animals or employ equipment until treatment has been absorbed, set or dried.

STORAGE AND DISPOSAL: Store airtight at room temperature. Nonrefillable container. Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard in trash.

SDS # CAD3212

Available exclusively from:

Contec, Inc. / Sporcidin Division

525 Locust Grove

Spartanburg, SC 29303 USA

1-864-503-8333

www.contecinc.com

2607-01/17



CONTEC® CITRIC ACID DISINFECTANT

BROAD SPECTRUM HOSPITAL DISINFECTANT

EPA Reg. No. 34810-35-71670

EPA Est. No. 34810-MO-1

ACTIVE INGREDIENTS

Citric Acid0.60%

OTHER INGREDIENTS99.40%

TOTAL 100.00%

For full list of ingredients, visit www.contecinc.com/saferchoice/citricaciddisinfectant

- Botanical
- Dye-free
- Fragrance-free
- Ready-to-use
- Germicidal
- Tuberculocidal
- Fungicidal
- Virucidal*

Keep Out of Reach of Children

1 U.S. Quart (32 fl. oz.)

Part No. CAD3212



CONTEC®



MADE IN USA